Gaussia luciferase

Catalog Number: 0308

Product Description

Gaussia luciferase is a luminescent protein with quantum yields higher than the luciferases from Renilla sp. and firefly. The enzyme catalyzes the oxidation decarboxylation of coelenterazine to give light (λ max = 485nm) and coelenteramide. ATP is not needed for activity. Gaussia luciferase is salt tolerant and retains it activity in the presence of non-ionic detergents, cholate, and deoxycholate. In addition, the enzyme is thermal tolerant and can retain activity at temperatures up to 60°C. The protein is also resistant to irreversible denaturation. Enzymatic activity can be recovered after treatment with 7M guanidine HCl or 8M urea in the presence of detergent.

The Avitag biotinylated version of the luciferase is no longer available.

Instructions for Use

- 1. The amount of luciferase used for a given assay should be empirically determined by titering the enzyme. The amount of luciferase that gives the highest signal: noise ratio should be selected.
- 2. Prepare a colenterazine stock by dissolving the powder in 1mL of acidified methanol (50µL concentrated HCl to 10mL of anhydrous methanol) to give a 200µM solution. Aliquot and store at -80°C for up to 4 weeks. A working solution of coelenterazine is prepared by diluting the stock solution to 10µM in 18Mohm water.

Important Notice: Coelenterazine is unstable in aqueous solutions. It is recommended that working solutions are prepared daily. Activity results will vary depending on concentration of Coelenterazine used.

3. For performing luciferase assays, dilute the enzyme in 50mM Tris-Cl, 0.5M NaCl, pH 7.8 to the desired RFU per assay. Add coelenterazine to 2.5µM and immediately measure the light production at 485nm with a 5-10 second integration.

Material Safety Data

FOR RESEARCH USE ONLY. NOT INTENDED OR APPROVED FOR HUMAN, DIAGNOSTICS OR VETERINARY USE. Do not ingest, swallow or inhale. Do not get in eyes, on skin, or on clothing. Wash thoroughly after handling. For complete safety information see full Material Safety Data Sheet.



Athena Enzyme Systems™ 1450 South Rolling Road Baltimore, MD 21227 USA

T (MD): 410-455-6319 T (USA): 888-892-8408 F: 410-455-1155 aesinfo@athenaes.com

Product Specifications

Protein	1mg (0308-1), 250µg (0308-2)
Long-term Storage	2 years at -80°C Aliquot to avoid repeated freezing and thawing
Short-tem Storage	1 month at 4°C
Formulation	In: 10 mM Na ₂ HPO ₄ , 140 mM NaCl, 2 mM KH ₂ PO ₄ , 3 mM KCl, 10% glycerol, pH 7.6
Molecular Weight	22kDa by SDS-PAGE





Figure 1. The graph on the top shows the emission spectrum for *Gaussia* luciferase (λ^{e}_{max} = 485nm). The lower graph shows the effects of salt concentration and pH on light emission of Gaussia luciferase.

References

^{1.} Verhaegen, M. and Christopoulos, T. K. 2002, Recombinant Gaussia Luciferase, Overexpression, Purification, and Analytical Application of a Bioluminescent Reporter for DNA Hybridization. Anal. Chem. 74: 4378-4385