

# Searching for **SILAC** Reagents?

## custom serum-free TISSUE CULTURE MEDIA



AthenaES™ offers custom cell culture media such as DMEM, RPMI, and EMEM for all research needs including SILAC (stable isotope labeling with amino acids in cell culture) experiments and animal-product-free applications. We can tailor the culture media formulation to your specifications, such as without specific amino acids like L-arginine and L-lysine.

Custom order volumes can be as little as 5L . If you would like to request a custom tissue culture media quote or place an order, call us at 888-892-8408, or submit an online form at [athenaes.com/custom\\_tissue.php](http://athenaes.com/custom_tissue.php).



**Athena Enzyme Systems™**  
1-888-892-8408  
1450 S. Rolling Rd. Baltimore, MD 21227  
[www.athenaes.com](http://www.athenaes.com)