

Product Data Sheet

Anti-YebF Antisera

Catalog Number: 0313

Product Description

Anti-YebF antisera is used for detecting the YebF protein in immunoassays. The antisera was prepared by injecting rabbits with purified YebF protein isolated from *E. coli* strain HMS174 harboring pAES40.

Instructions for Use

The antisera is supplied diluted 1:10 in 1x PBS pH 7.3. Store at 4°C. Do not store the antibody at the working concentration.

For use in Western blot, dilute to 1:100 to 1:500 in Tris or phosphate buffered saline pH 7.3 with 1% carrier protein (non-fat dry milk or BSA). The concentration of antibody needed for optimal results may vary and should be verified.

For other types of immunoassays, it is recommended that the antibody be titrated to determine the optimal concentration.



Athena Enzyme Systems™
1450 South Rolling Road
Baltimore, MD 21227
AthenaES®
a division of Athena Environmental Sciences, Inc.

T (MD): 410-455-6319
T (USA): 888-892-8408
F: 410-455-1155
USA
support@athenaes.com

Material Safety Data

FOR RESEARCH USE ONLY. NOT INTENDED OR APPROVED FOR HUMAN, DIAGNOSTICS OR VETERINARY USE. Do not ingest, swallow or inhale. Do not get in eyes, on skin, or on clothing. Wash thoroughly after handling. For complete safety information see full Material Safety Data Sheet.

Storage Conditions

Store product desiccated at -20°C.
Store reconstituted product at 4°C.